

# Effect of Next Generation Tobacco Products on Teeth Colouration

#1568



BRITISH AMERICAN TOBACCO

Annette Dalrymple<sup>1</sup>, Peter J Hall<sup>2</sup>, Thomas C Badrock<sup>2</sup>, Paul S Spradbery<sup>2</sup>,  
Matthew J Gornall<sup>2</sup>, James J Murphy<sup>1</sup> and Marianna D Gaca<sup>1</sup>

<sup>1</sup>British American Tobacco (Investments) Ltd, Southampton, SO15 8TL, UK


<sup>2</sup>Intertek Clinical Research Services, Hooton, Cheshire, CH66 7NZ, UK

Correspondence: Annette\_Dalrymple@bat.com

## Introduction

- Cigarette smoking is associated with changes in the oral cavity including tooth discolouration [1]
- Next generation tobacco and nicotine products (NGPs) such as e-cigarettes (e-cig) and tobacco heating products (THP) have rapidly increased in popularity over the last decade
- There is also a growing consensus that NGPs hold great potential for reducing the risk associated with cigarette smoking [2-3]
- However, the effects of e-cig and THP aerosols on the oral cavity or tooth discolouration are unknown

## Products



**Figure 1:** Products assessed for enamel discolouration *in vitro*. A) 3R4F Kentucky scientific reference cigarette (<https://ctcp.uky.edu/>), B) Glo™ commercial THP with Bright Tobacco Neostick™, C) Schematic of BAT Prototype e-cig, a: battery, b: consumable (containing e-liquid with 5 mg/mL nicotine and flavours), c: device body.

## Methods

### Total particulate matter (TPM) generation

- Products were attached to a linear smoke machine and smoke/aerosol captured on Cambridge filter pads (CFP)
- The puffing regimes used were:
  - 3R4F cigarette and Glo™: Health Canada Intense [4] - 55 mL puff volume, 2 sec puff duration, 30 sec puff interval. Filter ventilation was blocked on 3R4F products
  - Prototype e-cig: CORESTA recommended method No. 81 [5] - 55 mL puff volume, 3 sec puff duration, 30 sec puff interval
- Figure 2A details representative Cambridge filter pads with collected aerosol
- The smoke/aerosol captured on CFPs were eluted using DMSO (Figure 2B)



**Figure 2.** A) Representative CFP from 3R4F (10 puffs), Glo™ (48 puffs), e-cig (500 puffs) and blank. B) 3R4F TPM eluted in DMSO.

### Enamel block preparation

- Enamel blocks (6.5 x6.5 mm) were cut from bovine incisors and polished using 400-grit silicon carbide paper to produce a surface finish similar to human enamel ( $R_a \approx 0.5 \mu\text{m}$ )
- Blocks were sonicated in deionised water for 2 minutes and then stored in phosphate-buffered saline (PBS)
- Blocks were incubated with sterile human saliva for 3 hours at 35°C to form a pellicle layer

### Enamel block exposure to TPM

- 10 enamel blocks were prepared per TPM and DMSO control
- Samples were incubated at 35°C for 14 days

### Colour measurement


- Blocks were periodically removed from the test article and hydrated in PBS for 60 minutes
- Colour readings ( $L^*$ ,  $a^*$ ,  $b^*$ ) were measured using a Konica Minolta CM-700d at 0, 1, 2, 3, 4, 7, 8, 9, 10, 11 and 14 days
- Discolouration was determined as change in colour ( $\Delta E$ ) between baseline and exposure day ( $\Delta E = \sqrt{((\Delta L^*)^2 + (\Delta a^*)^2 + (\Delta b^*)^2)}$ )
- A one-way ANOVA test was used to test for a difference in  $\Delta E$  means at Day 14 between the four test articles. The Tukey procedure was used for pairwise comparisons

## Results

- Discolouration was apparent in as little as 1 day exposure to 3R4F TPM and continued to increase until day 14 (Figures 3 and 4)
- The mean colour change induced by 3R4F TPM at Day 14 was statistically significantly ( $p < 0.0001$ ) higher than other treatments
- Mean  $\Delta E$  values were  $29.39 \pm 3.57$ ,  $10.67 \pm 2.64$ ,  $10.46 \pm 2.30$  and  $12.57 \pm 1.97$  for 3R4F, e-cig, THP and DMSO respectively



Annette is a Claims Implementation Senior Scientist at British American Tobacco (BAT) based in Southampton UK. She has a PhD in reproductive biology from Edinburgh University and went on to postdoctoral research at King's College London for 6 years before joining BAT R&D in 2007. Annette has been instrumental in the development and application of new laboratory methods to investigate next generation nicotine and tobacco products (NGPs). She is an active member of the tobacco scientific community, a member of the Society of Toxicology and has published many papers on the development and application of methods for assessing the toxicity of tobacco products and more recently e-cigarettes. Please contact: [annette\\_dalrymple@bat.com](mailto:annette_dalrymple@bat.com)



**Figure 3.** Mean  $\Delta E$  values of enamel blocks at days 0-14 for 3R4F cigarette, BAT prototype e-cig and Glo™ TPM or DMSO control.



**Figure 4.** Enamel block discolouration after 14 days exposure to TPM generated from 3R4F cigarettes, Glo™ or BAT prototype e-cig and also DMSO control.

## Conclusions

- Cigarette smoke (3R4F) significantly increased the level of discolouration on bovine enamel teeth samples
- E-cig and THP exposure produced little discolouration and were comparable to the 'untreated' DMSO control group
- Further studies are required to assess the longer-term impact of e-cig and THP aerosols on tooth discolouration

## References

1. M.N. Alkhatib, R.D. Holt, R. Bedi. Smoking and tooth discolouration: findings from a national cross-sectional study. BMC Public Health. 2005 Mar 24;5:27
2. A. McNeill, L.S. Brose, R. Calder, S.C. Hitchman, P. Hajek, H. McRobbie, E-cigarettes: an evidence update. A report commissioned by Public Health England. [https://www.gov.uk/government/uploads/system/uploads/attachment\\_data/file/457102/E-cigarettes\\_an\\_evidence\\_update\\_A\\_report\\_commissioned\\_by\\_Public\\_Health\\_England\\_FINAL.pdf](https://www.gov.uk/government/uploads/system/uploads/attachment_data/file/457102/E-cigarettes_an_evidence_update_A_report_commissioned_by_Public_Health_England_FINAL.pdf), 2015 (accessed 18.07.17).
3. Royal College of Physicians, Nicotine without smoke: tobacco harm reduction. <https://www.rcplondon.ac.uk/projects/outputs/nicotine-without-smoke-tobacco-harm-reduction-0>, 2016 (accessed 18.07.17).
4. Health Canada Official Method T-115. Determination of "Tar", Nicotine and Carbon Monoxide in Mainstream Tobacco Smoke. Prepared by the Department of Health, Dec. 31, 1999.
5. CORESTA, No. 81 - routine analytical machine for e-cigarette aerosol generation and collection - Definitions and standard conditions. [https://www.coresta.org/sites/default/files/technical\\_documents/main/CRM\\_81.pdf](https://www.coresta.org/sites/default/files/technical_documents/main/CRM_81.pdf), 2015 (accessed 18.07.17)

[www.bat-science.com](http://www.bat-science.com)



@BAT\_Sci

