## Platelet and Monocyte Binding to Human Aortic Endothelial Cells in well plate microfluidic devices after TNFα stimulation

Om Makwana<sup>1</sup>, Natalia Cockcroft<sup>2</sup>, Hannah Flockton<sup>3</sup>, Gary Watters<sup>3</sup>, Rizwan Nisar<sup>3</sup>, Gina Smith<sup>3</sup>, Wanda Fields<sup>1</sup> and Betsy Bombick<sup>1</sup>

<sup>1</sup>RAI Services Company, Winston-Salem, NC USA

<sup>2</sup>BAT GR&D Southampton, UK

<sup>3</sup>Covance Laboratories Ltd, Genetic and Molecular Toxicology, Harrogate, UK

# Abstract #3347

#### Abstract

Well plate microfluidic devices have been demonstrated to provide a physiologically relevant environment to model the human vasculature. Here we utilized the BioFlux™ microfluidic technology in conjunction with human aortic endothelial cells (HAECs) to observe the development of TNFα-induced platelet and monocyte adhesion to endothelial cells (i.e., thrombosis). HAEC monolayers were grown in microfluidic channels and were then exposed to a range of TNFα concentrations (5-100 ng/mL) for 5 hours under gravity flow. This was followed by perfusion through the microfluidic channel with whole blood (preincubated TNFα for 1 hour) at 5 dyn/cm<sup>2</sup> for 5 minutes or THP-1 cells (activated with TPA) at 1 dyn/cm<sup>2</sup> for 1 minute followed by 0.5 dyn/cm<sup>2</sup> for 10 minutes. Whole blood and THP-1 cells were pre-labeled with Calcein AM to aid quantification of cell binding (percent channel area) with BioFlux™ Montage software. We demonstrate that both platelets and monocytes (THP-1) are able to bind to HAECs in microfluidic devices after exposure to TNFα and furthermore we have shown that this trend appears to increase in a dose-responsive manner before plateauing in both whole blood and THP-1 experimentation. We believe that these microfluidic devices have the potential to evaluate pro-inflammatory compounds and mixtures that may be relevant to cardiovascular disease such as atherosclerosis.

#### Material and Methods

HAECs were maintained in VascuLife® VEGF Endothelial Cell Culture Medium supplemented with VEGF LifeFactor Kit (Lifeline Cell Technology). Human THP-1 monocytic leukemia (THP-1) cells were supplied by ATCC and maintained in the recommended RPMI-1640 medium. Fresh whole blood samples were collected from fasting donors using 0.32% sodium citrate BD vacutainers™ and were used within 3 hours of collection.

Methodology for BioFlux™ system and plate set up were adapted from Conant et al., 2009 and 2011. HAEC monolayers were exposed to TNFα under gravity flow (<0.1 dyn/cm²) for 5 or 24 hours as required. Whole blood was labelled with 4 μM Calcein AM (0.5% v/v). TNFα was incubated with whole blood (10% v/v) for 1 hour on a slow moving rocking platform. The whole blood was perfused through channels at 5 dyn/cm<sup>2</sup> for 5 mins and followed by a wash with HBSS containing 10 µg/mL Hoechst 33342 (H33342) for 1 min.

To activate THP-1 cells, 12-O-tetra-decanoylphorbol-13-acetate (TPA) in DMSO was added to cell samples (<0.1% v/v) at 1 ng/mL and transferred to low adhesive and adhesive flasks for 24 hours to determine successful activation. THP-1 cells were not pre-incubated with TNFα, but were labelled with 4 μM Calcein AM (0.5% v/v). For 5-hour experiments, THP-1 cells were perfused through channels at 1 dyn/cm<sup>2</sup> for 1 minute followed by 0.5 dyn/cm<sup>2</sup> for 10 min. Channels were washed with HBSS containing H33342.

Images were acquired using the Bioflux™ Montage software from multiple fields of view per channel under the brightfield, FITC (green; Calcien AM) and DAPI (blue; H33342) filters. Platelet/THP-1 adhesion was calculated using Montage software based on pixel area. Oneway ANOVA with multiple comparisons \*\*\* p<0.001.

For 24-hour TNFα experimentation, un-labelled, activated THP-1 cells were perfused at 2 dyn/cm<sup>2</sup> for 30 seconds followed by 0.5 dyn/cm<sup>2</sup> for 1 min. THP-1 cells were left for a static adhesion period followed by a wash with HBSS containing H33342 at 2, 5 and 6 dyn/cm<sup>2</sup> for 1 min. The perfusion, adhesion period and wash cycle was performed four times with increasing adhesion periods (10, 20, 30 and 40 mins). THP-1 numbers were quantified by manual counting of phase contrast images from three fields of view. Two-way ANOVA with multiple comparisons p<0.001.

#### Summary and Conclusions

- > Both platelets and monocytes are able to bind to HAECs in microfluidic devices after TNFα exposure for 5 hours (Figures 1 & 2).
- > Binding of monocytes to HAECs in microfluidic devices is increased after TNFα exposure for 24 hours (Figure 3).




Figure 1A. Whole Blood adhesion to HAECs following TNFα stimulation (1-hour Blood + 5-hour HAECs) \*\*\* p<0.001




Figure 1B. Representative images of Whole Blood (FITC) adhesion to HAECs (DAPI) following 5-hour TNFα incubation




Figure 3. THP-1 Adhesion to HAECs following repeat THP-1 addition (24 hour TNFα)




Figure 2A. THP-1 adhesion to HAECs following TNFα stimulation for 5-hours \*\*\* p<0.001




Figure 2B. Representative images of THP-1 (FITC) Adhesion to HAECs (DAPI) following 5-hour TNFα incubation

### References

- Conant CG, Schwartz MA, Nevill T, Ionescu-Zanetti C (2009). Platelet adhesion and aggregation under flow using microfluidic flow cells. J Vis Exp. Oct 27;(32).
- Conant CG, Nevill JT, Zhou Z, Dong J-F, Schwartz MA, Ionescu-Zanetti C (2011). Using well-plate microfluidic devices to conduct shear-based thrombosis assays. J Lab Autom. Apr;16(2):148-52.
- Cockcroft NY, Oke O, Cunningham F, Bishop E, Fearon IM, Zantl R, et al. (2009). An in vitro perfusion system to examine the responses of endothelial cells to simulated flow and inflammatory stimulation. Altern Lab Anim. Dec;37(6):657–69.
- Tsai M, Kita A, Leach J, Rounsevell R, Huang JN, Moake J, et al. (2012). In vitro modeling of the microvascular occlusion and thrombosis that occur in hematologic diseases using microfluidic technology. J Clin Invest Jan;122(1):408-18.